

MicroFlask

***Miniaturize without compromise
Simple and low cost bioreactor scale-down to micro-liter volumes***



Introduction


The MicroFlask is a cultivation system that turns microtiter plates into reliable and easy to operate cultivation systems offering reproducible cultivation results.

Miniaturization of cultures to standard 24 or 96 well microtiter plate format is an attractive format because of the common availability of:

- Low-cost microtiter plates in various shapes and materials
- Robotic equipment for filling wells with growth medium
- Centrifugation (pelleting cells)
- Equipment for UV-VIS or fluorescence measurement, etc.


Main features of the MicroFlask:

- ✓ 24 (or 96) single use and low cost bioreactors per titer plate
- ✓ Stackable design allows up to 4 plates in one clamp. This results in up to 384 bioreactors on 85 x 125 mm
- ✓ Oxygen transfer rates are similar to those in Erlenmeyer flask (40 mmol O₂/l/min) and can be achieved with standard orbital shaking equipment
- ✓ No well-to-well cross-contamination risk
- ✓ Low (and uniform) evaporation rates
- ✓ Cell densities of 10 g dry wt/l or more are achievable
- ✓ Well-to-well reproducibility
- ✓ Sandwich covers converting each well of 24 and 96 well micro-titer plates (both deep-well and low-well plates) into individual "micro-reactors"
- ✓ A spring-loaded replicator allows the simultaneous and reproducible inoculation of MicroFlask cultures and sampling from 96 frozen glycerol stocks without thawing the main culture. The stocks can remain viable for years.



Examples of MicroFlask application areas:

- Discovery of secondary metabolites or new enzyme activities (heterogeneous culture collections)
- Screening and distribution of mutant and construct libraries e.g. in E. coli or yeasts
- Metabolic flux studies and high-throughput screening for high-activity mutants
- Comparative studies, e.g. clinical isolates
- Growth medium optimization for cell lines or production strains




Oxygen supply

Classic bioreactors for microbial cultures require an air supply of up to two working volumes per minute as an average. The sandwich covers are designed to have a headspace refreshment rate of 1-2 culture volumes per minute. This headspace refreshment rate secures oxygen concentrations in the headspace of at least 18% (v/v), even at an oxygen consumption rate of 40 mmol O₂/L/ minute. At the same time this limits the evaporation to about 2% of the culture volume per day (at an ambient humidity of 50% and 30°C).

Sandwich cover

- A stainless steel lid for rigid support
- B micro fiber layer
- C extruded teflon or ePTFE (0.2 µm pores) layer for sterility
- D stainless steel foil with pinholes to further limit of evaporation (only in covers for 96-lowwell plates)
- F soft silicone layer (to seal the “mini-reactors”)




OTR Specifications MicroFlask-bioreactors


Type of microtiter plate	Well volume dimension in horizontal plane	Culture volume	Orbital shaking frequency	Shaking amplitude	O ₂ -transfer rate (30°C, air, 1 bar)	Headspace refreshment rate	Evaporation rate per well (at 30°C)	Mixing pattern at 300 rpm	
								ampl. 25 mm	ampl. 30mm
24-square deep-well polypropylene, 17x17 mm, depth 40 mm		2500 µl	300 rpm	50 mm	51 mmol O ₂ / l	2.5 ml / min	50% humidity: 50 µl H ₂ O per day		
		2500 µl	300 rpm	25 mm	39 mmol O ₂ / l	(1 VVM)			
		2500 µl	220 rpm	50 mm	35 mmol O ₂ / l				
		4000 µl	300 rpm	50 mm	24 mmol O ₂ / l	2.5 ml / min	75% humidity: 25 µl H ₂ O per day	2500 µl	2500 µl
4000 µl	220 rpm	25 mm	24 mmol O ₂ / l	(0.6 VVM)					
24-round low-well polystyrene, Ø 16 mm, depth 18 mm		750 µl	300 rpm	50 mm	40 mmol O ₂ / l	1.1 ml / min	50% humidity: 30 µl H ₂ O per day		
		750 µl	300 rpm	25 mm	25 mmol O ₂ / l	(1.4 VVM)			
		1000 µl	300 rpm	50 mm	30 mmol O ₂ / l	1.1 ml / min	75% humidity: 15 µl H ₂ O per day	1000 µl	1000 µl
		1000 µl	300 rpm	25 mm	19 mmol O ₂ / l	(1.1 VVM)			
96-square deep-well polypropylene, 8x8 mm, depth 40 mm		500 µl	300 rpm	50 mm	38 mmol O ₂ / l	1 ml / min	50% humidity: 22 µl H ₂ O per day		
		500 µl	300 rpm	25 mm	12 mmol O ₂ / l	(2 VVM)			
		750 µl	300 rpm	50 mm	24 mmol O ₂ / l	1 ml / min	50% humidity: 22 µl H ₂ O per day	750 µl	750 µl
		750 µl	300 rpm	25 mm	7 mmol O ₂ / l	(1.3 VVM)			
		1000 µl	300 rpm	50 mm	18 mmol O ₂ / l	1 ml / min	75% humidity: 11 µl H ₂ O per day	750 µl	750 µl
		1000 µl	300 rpm	25 mm	3 mmol O ₂ / l	(1 VVM)			
96-round low-well polypropylene, Ø 6.5 mm, depth 10 mm		100 µl	300 rpm	50 mm	39 mmol O ₂ / l	250 µl / min	50% humidity: 6 µl H ₂ O per day		
		100 µl	300 rpm	25 mm	20 mmol O ₂ / l	(2.5 VVM)			
		150 µl	300 rpm	50 mm	32 mmol O ₂ / l	250 µl / min	50% humidity: 6 µl H ₂ O per day	150 µl	150 µl
		150 µl	300 rpm	25 mm	16 mmol O ₂ / l	(1.7 VVM)			
		200 µl	220 rpm	50 mm	12 mmol O ₂ / l	250 µl / min	75% humidity: 3 µl H ₂ O per day	150 µl	150 µl
		200 µl	300 rpm	25 mm	12 mmol O ₂ / l	(1.3 VVM)			

MicroFlask


Sandwich covers for 24-square deep-well plates allow a culture volume of 2.5 - 4 ml, and are often used for (pellet-forming) streptomycetes and filamentous fungi, e.g. for medium optimization studies or product discovery.




Sandwich covers for 96-square deep-well plates allow a culture volume of 0.5 - 1 ml and are most often used for bacterial and yeast cultures.




Sandwich covers for 24-low-well plates allow culture volumes of 0.5 - 1 ml and are especially suitable for eukaryotic cell cultures.



Sandwich covers for 96-low-well plates allow a culture volume of 0.1 - 0.2 ml and are most often used for bacterial mutant and construct libraries.



The cover clamps for deep-well plates
The cover clamp can be mounted on orbital shaking platforms from a range of suppliers.



The cover clamps for low-well plates
allow the stacking of up to three 24-well plates or four 96-well plates (hence a maximum of 12 or 16 plates per full clamp). The cover clamp can be mounted on **orbital** shaking platforms from a range of suppliers.

Literature list • Duetz, W.A., Rüedi, L., Hermann, R., O'Connor, K., Büchs, J., and Witholt, B. (2000). Methods for intense aeration, growth, storage, and replication of bacterial strains in microtiter plates. *Applied and Environmental Microbiology* 66:2641-2646. • Minas, W., Bailey, J.E., Duetz, W.A. (2000). Streptomycetes in microcultures: growth, production of secondary metabolites, and storage and retrieval in the 96-well format. *Antonie van Leeuwenhoek*, 78 (3-4): 297-305 • Duetz, W.A. and Witholt, B. (2001). Effectiveness of orbital shaking for the aeration of suspended bacterial cultures in square deepwell microtiter plates. *Biochemical Engineering Journal*, 7: 113-115 • Duetz, W.A., Minas, W., Kuhner, M., Clerval, R., Fjällman, A.H.M., Witholt, B. (2001). Miniaturized microbial growth systems in screening. *Bioworld* 2: 8-10 • Duetz, W.A., Witholt, B. (2004). Oxygen transfer by orbital shaking of square vessels and deepwell microtiter plates of various dimensions. *Biochemical Engineering Journal* 17: 181-185



Applikon Biotechnology BV De Brauwweg 13, P.O. Box 149, 3100 AC Schiedam, The Netherlands.

Phone: +31 10 298 35 66. Fax: +31 10 437 96 48. E-mail: biotech@applikon.com

www.applikon.com