

TECRA®

MICROBIOLOGICAL ENRICHMENT MEDIA & SUPPLEMENTS

A range of high quality Enrichment Media, including proprietary formulations and Selective Supplements, specially developed to maximise the performance of TECRA's diagnostic test kits.

Why use TECRA media with TECRA kits?

Proven Performance

- ◆ TECRA media have been **fully validated** for use with TECRA's diagnostic kits.
- ◆ Proprietary media formulations have been developed and proven to optimize performance.
- ◆ This gives you confidence that TECRA media will give superior results.

Convenience

- ◆ TECRA media are formulated to minimise preparation time and thereby save you money.
- ◆ TECRA media can be conveniently supplied together with TECRA kits.

Customer and Product Support

- ◆ Comprehensive product training courses on the use of TECRA products are available to ensure that laboratory and quality assurance staff are competent in the use of TECRA methods.
- ◆ TECRA has a commitment to excellent service so you can rely on extensive technical and product support either directly from TECRA or through TECRA's global network of distributors.

Quality & Consistency

- ◆ TECRA media are manufactured from premium quality raw materials and to the highest quality standards.
- ◆ Each batch of TECRA media is tested to ensure it meets stringent quality assurance specifications and batch to batch consistency.
- ◆ All TECRA media products carry a quality certificate. This is your guarantee of quality.
- ◆ You can expect consistent performance time after time by using TECRA media which have been previously validated for use with your chosen test method.



TECRA offers everything you need to achieve reliable results:

- ◆ Validated, quality test kits
- ◆ Media formulated for use with the kits
- ◆ Proven and simple to use methods to optimize results
- ◆ Technical support and training

TECRA Media & Media Supplements

MEDIA	Product Code	Unit
TECRA Buffered Listeria Enrichment Broth	BLEMED500	500g
Buffered Peptone Water	BPWMED500	500g
Buffered Peptone Water	BPWMED2000	2000g
TECRA Campylobacter Enrichment Broth Base (selective agents supplied separately)	TCE(B)MED500	500g
Fraser Broth Base (selective agent supplied separately)	FRB(B)MED500	500g
TECRA Lactose Broth	TLBMED500	500g
TECRA Listeria Enrichment Broth (selective agents included)	TLEMED500	500g
TECRA Listeria Enrichment Broth Base (selective agents supplied separately)	TLE(B)MED500	500g
M Broth	TMBMED500	500g
Modified Buffered Peptone Water	MBPMED500	500g
Modified Buffered Peptone Water	MBPMED2000	2000g
Modified EC Broth	MECMED500	500g
Modified Lethen Broth	MLBMED500	500g
TECRA Pseudomonas Enrichment Broth Base (selective agents supplied separately)	TPE(B)MED500	500g
Rappaport-Vassiliadis R10 Broth	RVRMED500	500g
TECRA Staphylococcus Growth Medium	TSGMED500	500g
Tetrathionate Broth Base (selective agents supplied separately)	TTB(B)MED500	500g
Tryptone Soy Broth	TSBMED500	500g
University of Vermont Medium	UVMMED500	500g
MEDIA SUPPLEMENTS	Product Code	Unit
TECRA Campylobacter Selective Supplement (for use with TECRA Campylobacter Enrichment Broth Base)	TCESUP010	10 vials
Imbentin Supplement	IMBSUP50	50mL
Imbentin Supplement	IMBSUP1000	1L
Imbentin Supplement	IMBSUP2000	2L
TECRA Listeria Selective Supplement (for use with TECRA Listeria Enrichment Broth Base)	TLESUP010	10 vials
TECRA Pseudomonas Selective Supplement (for use with TECRA Pseudomonas Enrichment Broth Base)	TPESUP010	10 vials
TECRA Salmonella Supplement	SALSUP250	250mL
TECRA Salmonella Supplement	SALSUP500	500mL

* For more information on TECRA media and supplements, refer to the TECRA methods manuals, the TECRA product list, visit the TECRA website or contact TECRA.

Your Partner in Rapid Testing



www.tecra.net

TECRA International Pty Ltd
 13 Rodborough Road, Frenchs Forest NSW 2086, AUSTRALIA
 Telephone: +61 2 8977 3000 Facsimile: +61 2 9453 3422
 E-mail: enquiries@tecra.net

Distributed by:

TECRA[®]
Excellence in Rapid Testing